

SARS-CoV-2 Antigen Test

Rapid SARS-CoV-2 Antigen Test (Nasopharyngeal specimen)



Step 1: Prepare all items for testing.

Step 2: Confirm product name and package.

Step 3: Check the expiry date. **Step 4:** Write an identifier on the cassette.



Step 5: Take a new disposable swab and unwrap it.

Step 6: Insert the swab through one nostril parallel to the palate (not upwards) until resistance is encountered.

Step 7: Gently rub and roll the swab over surface of the nasopharynx for 10 times to absorb secretions. Withdraw the swab from the nasal cavity slowly.



Step 8: Add 7 drops (~300µl) of buffer into the test tube.

Step 9: Insert the tip of the swab into the solution, and rotate it close to tube inner wall about 10 times.

Step 10: Squeeze the tip of the swab along the tube inner wall, then remove and discard the swab.

Step 11: Cover the test tube with the cap.



Step 12: Add 3 drops (~100µl) of treated specimen into "S" well.

Step 13: Interpret the result between 15 to 20 minutes.

